BLUE

**AESTHETIC**

POLYMER/MONOMER

**INSTRUCTIONS FOR USE**

---

GEBRAUCHSINFORMATION

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---

Dear customer,

The product that you have purchased qualifies as a medical device under European Directive 93/42/EEC. Please ensure that you file the lot number and product name in your administrative system on receipt of the product. For each piece of work that you produce, please make a note of all the materials used along with all the appropriate lot numbers. For your own protection as well as for the protection of your patients and the environment, please also follow the important guidelines in the following instructions for use.

#### **Product description**

AESTHETIC BLUE cold-curing denture base material has excellent flow and modelling properties which make it easy to use with either the pouring or packing technique. AESTHETIC BLUE cold-curing denture base material is reliable even when there are two or more multiple denture base segments. The dental technician has different opaque shades, transparent pink and clear to choose from. The variable mix proportions of the powder and liquid within normal limits allow the consistency and flexibility in handling to be determined on an individual basis.

#### **Composition / Powder**

Polymethyl methacrylate, softening agent, benzol peroxide, catalyst, pigments

EN

---

### **Composition / Liquid**

Methyl methacrylate, dimethacrylate, catalyst

### **Storage instructions**

Store the material in a cool, dark, well-ventilated place.

Storage temperature: 2–28 °C / 36–82 °F. Do not use the material after the expiry date.

Keep out of the reach of children.

### **Indication**

- Full dentures
- Partial dentures
- Combination dentures
- Dental implants
- Repairs

### **Contraindication**

Avoid direct contact with unpolymerised material within the oral cavity. If a patient is known to be allergic to any of the ingredients in AESTHETIC BLUE, the material must not be used.

### **Side effects**

No systemic side effects have been reported to date. In individual cases, local allergic reactions to PMMA/MMA-based denture base materials have been reported.

### **Danger warnings**

- Monomer contains methyl methacrylate (MMA)
- MMA is an irritant and easily flammable (flash point: + 10 °C / 50 °F)
- MMA and its vapours are irritating to the eyes, skin and respiratory system
- May cause sensitisation upon skin contact
- Do not inhale vapours
- Keep away from sources of ignition – no smoking

- Do not empty into drains
- Avoid contact of the skin with monomer and uncured material. Many commercial gloves, e.g. those made of latex or vinyl, are not monomer-resistant and therefore do not provide protection against the sensitising effect of methacrylates.
- Wear a mask when grinding and use a suction removal system.
- Always use a spatula when handling the mixture.
- Take precautionary measures against static discharges.

## Warning

This material has been developed solely for use in dentistry and must be handled strictly in accordance with the instructions for use. The manufacturer cannot be held liable for any damages arising as a result of failure to observe the instructions for use or the stated area of application. The same applies in the event that the product is mixed or processed with other manufacturers' products. The user shall be solely responsible for testing the material with respect to its suitability prior to use for any purpose other than those explicitly stated in the instructions.

**Safety data sheets can be found on our website at [www.candulor.com](http://www.candulor.com).**

## INSTRUCTIONS FOR USE

---

### Important

- The precast and model must be well soaked in water.
- Apply two coats of ISO-K separating liquid and allow to dry thoroughly.
- Roughen the teeth and wet with monomer to ensure proper bonding with the denture base.
- The length of the pouring and modelling phases will depend on the chosen mixing ratio and room temperature.
- Mix the powder and liquid mixture thoroughly using a hard spatula (approx. 20–30 sec.).
- Avoid cooling quickly in cold water (causes stress cracks)

## POURING TECHNIQUE

---

### 1. Preparation


After boiling out the water-saturated plaster, isolate the surfaces with two applications of ISO-K separating liquid and allow to dry thoroughly. Roughen the teeth well using a dental drill and wet with monomer in order to ensure proper bonding with the denture base.

### 2. Dosage

Ideal mixing ratio for one denture:

15 g polymer : 10 ml monomer

**PolyMaster:** 30 g polymer : 15 ml monomer

The dosage system ensures an ideal mixing ratio and minimum polymerisation shrinkage of AESTHETIC BLUE. The scale on the monomer cylinder is in millilitres, on the polymer beaker in grams.

### Dosage as required within normal limits

Pour the required quantity of monomer into the plastic cup and add an appropriate amount of polymer to obtain a viscous consistency. Too high a monomer content will alter the shade, consistency and handling properties of the material. This will also inhibit polymerisation, possibly leading to irritation of the mucous membranes.

### 3. Mixing

Mix the powder and liquid thoroughly in the mixing ratio or as required and blend thoroughly for approx. 20–30 seconds using the spatula. Leave to stand for 15 seconds so that any air bubbles can escape.

**Mix the powder and liquid thoroughly for approx. 20–30 seconds!**

### 4. Pouring phase

The pouring phase lasts approximately 2,5–3 minutes at room temperature (18–25 °C / 64–77 °F). Pour the material into the precast within this time span.

### 5. Modelling phase

After a transition phase of roughly 4 minutes the material will assume a firmer consistency and can be modelled for a period of 3 minutes. The higher the room temperature, the shorter the pouring and modelling phases!

**The synthetic material must be placed in the pressure vessel within no more than 8–9 minutes.**

## 6. Polymerisation

Polymerisation takes place in the pressure vessel for 15 minutes at a temperature of 40 °C / 104 °F and a pressure of 2 bar (check).

## 7. Finishing

After removing the precast, finish in the usual manner.


## PACKING TECHNIQUE

### 1. Preparation

After boiling out the water-saturated plaster, isolate the surfaces with two applications of ISO-K separating liquid and allow to dry thoroughly. Roughen the teeth well and wet with monomer to ensure proper bonding with the denture base.

### 2. Dosage

Ideal mixing ratio for one denture:  
20,5 g polymer : 10 ml monomer

The dosage system ensures an ideal mixing ratio and minimum polymerisation shrinkage of AESTHETIC BLUE. The scale on the monomer cylinder is in millilitres, on the polymer beaker in grams.

Too high a monomer content will alter the shade, consistency and handling of the material. This may also inhibit polymerisation, possibly leading to irritation of the mucous membranes and to porosity.

### **3. Mixing**

Mix the powder and liquid thoroughly in the mixing ratio for approx. 20–30 seconds using the spatula. Put the lid on the mixing beaker and leave to stand for roughly 3 minutes (at a room temperature of 23 °C / 73 °F). The resulting mixture is then workable for a period of approximately 2 minutes. The higher the room temperature, the shorter the working and setting times.

### **4. Pressing**

Place an ample amount of the mixture in one half of the flask which must be at body temperature and which you have previously wetted and isolated with ISO-K. Close the flask carefully, place it under the press and apply a pressure of 80 bar, or place it in the clamp.

## 5. Polymerisation

Polymerisation takes place in the clamp or press under constant pressure over a period of 30 minutes (at a room temperature of 23 °C / 73 °F).

## 6. Removal and finishing

Open the flask and remove the plaster. Check the occlusion of the denture and finish in the usual manner.

## Repairs and corrections using AESTHETIC BLUE

Repairs and corrections to AESTHETIC BLUE may be carried out using AESTHETIC BLUE and the pouring technique. The contact surfaces must be well roughened and wetted with monomer.

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🏢 Candulor AG