# AESTHETIC

COLOR SET EASY

**INSTRUCTIONS FOR USE**

---

GEBRAUCHSINFORMATION

## INHALTSVERZEICHNIS / CONTENTS

---

		SEITE / PAGE
DE	Gebrauchsinformation	4
EN	Instructions for use	24
FR	Instructions de mise en oeuvre	42
IT	Istruzioni d'uso	62
ES	Instrucciones de uso	82
NL	Gebruiksaanwijzing	102
RU	Инструкция по применению	122

Manufactured for Candulor  
**Rx only – For dental use only!**  
Complies with/entspricht: ISO 20795-1; EN ISO 20795-1

DE

---

EN

---

FR

---

IT

---

ES

---

NL

---

RU

---

## AESTHETIC COLOR SET EASY

### AESTHETIC RED/AESTHETIC BLUE

---

#### Inhalt

- 4 x 100 g AESTHETIC RED bzw. AESTHETIC BLUE Polymer
- 1 x 150 ml AESTHETIC RED bzw. AESTHETIC BLUE Monomer
- 1 x 150 ml AESTHETIC LT Modelling Monomer (RED bzw. BLUE)
- 1 x Fibrillen
- 8 x 15 g AESTHETIC Intensive Colors
- 1 x 15 g AESTHETIC Intensive Opaque

DE

---

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---

#### Contents

- 4 x 100 g AESTHETIC RED resp. AESTHETIC BLUE Polymer
- 1 x 150 ml AESTHETIC RED resp. AESTHETIC BLUE Monomer
- 1 x 150 ml AESTHETIC LT Modelling Monomer (RED resp. BLUE)
- 1 x Fibers
- 8 x 15 g AESTHETIC Intensive Colors
- 1 x 15 g AESTHETIC Intensive Opaque

- 1 x Shade guide
- 1 x Modelling Instrument
- 1 x Glass pipette
- 1 x Ceramic Mixing Plate
- 1 x Instructions for use AESTHETIC COLOR SET EASY
- 1 x Instructions for use AESTHETIC RED resp. AESTHETIC BLUE
- 1 x Various auxiliary materials

## AESTHETIC COLOR SET EASY MINI

### AESTHETIC RED/AESTHETIC BLUE

---

#### Contents

- 3 x 100 g AESTHETIC RED resp. AESTHETIC BLUE Polymer
- 1 x 150 ml AESTHETIC LT Modelling Monomer (RED resp. BLUE)
- 1 x Shade guide
- 1 x Glass pipette
- 1 x Instructions for use AESTHETIC COLOR SET EASY
- 1 x Instructions for use AESTHETIC RED resp. AESTHETIC BLUE
- 1 x Various auxiliary materials

## AESTHETIC INTENSIVE COLORS

---

#### Contents

- 1 x Snow White (01)
- 1 x Ocean Blue (02)
- 1 x True Pink (03)
- 1 x Purpur Red (04)
- 1 x Sunny Orange (05)
- 1 x Mystic Braun (06)
- 1 x Dark Shadow (07)
- 1 x Bright Light (08)
- 1 x Opaque

## Composition

### **AESTHETIC RED Monomer / AESTHETIC LT Modelling Monomer**

- > 85 % Methyl methacrylate
- > 5 % Ethyleneglycol-dimethacrylate

### **AESTHETIC BLUE Monomer / AESTHETIC LT Modelling Monomer**

- > 90 % Methyl methacrylate
- < 5 % Butanediol-dimethacrylate

### **AESTHETIC Intensive Opaque / AESTHETIC Intensive Colors**

- > 99 % Polymethyl methacrylates
- < 0,3 % Pigments (ferrous oxide)
- 0,5–1,5 % Benzoyl peroxides

### **AESTHETIC Intensive Color white**

- > 98 % Polymethyl methacrylates
- < 1 % Pigments (titanium dioxide)

### **Fibers**

Do not contain hazardous substances in the meaning of EU Directive 67/548/EEC.

Dear customer,

The product that you have purchased qualifies as a medical device under European Directive 93/42/EEC. Please ensure that you file the lot number and product name in your administrative system on receipt of the product. For each piece of work that you produce, please make a note of all the materials used along with all the appropriate lot numbers. For your own protection as well as for the protection of your patients and the environment, please also follow the important guidelines in the following instructions for use.

## PRODUCT DESCRIPTION

---

### **AESTHETIC LT Modelling Monomer**

The respective AESTHETIC LT Modelling Monomer for AESTHETIC BLUE or AESTHETIC RED is used to characterise the denture. With AESTHETIC LT Modelling Monomer we obtain a longer modelling phase which also prevents the acrylic from drying out.

### **Polymer**

Ready-mixed colors in three different shades enable you to produce individual coloring of the gingiva simply and above all quickly.

EN

---

### **AESTHETIC Intensive Opaque / AESTHETIC Intensive Colors**

With AESTHETIC Intensive Opaque a higher degree of opacity is achieved. This can reduce the visibility of metal parts particularly in the case of partial or implant dentures. Another indication is the covering of thin layers of acrylic in the area of the necks of the teeth. (Ratio 1 : 20)

AESTHETIC Intensive Colors can be used for individual shading of the denture base materials. The coloring can be reproduced with the aid of the sketch.

### **Fibres**

The red fibres are directly added to the monomer or can be mixed with the acrylic. Fibres (veins) which have been mixed in with the monomer are added after forming the layers of the alveolar mucosa.

### **Preparation**

In view of the difference in dough time between the AESTHETIC denture base materials (AESTHETIC RED 8–10 min./ AESTHETIC BLUE 3 min.) the time for mixing is:

- **before** characterisation in the case of heat-curing denture base material
- **after** characterisation in the case of cold-curing denture base material

For this purpose we use the usual AESTHETIC monomers for cold-curing or heat-curing base material.

## INSTRUCTIONS FOR USE

---

1. To start off, we mix the basic shade (polymer) 53/55/57 with AESTHETIC LT Modelling Monomer. As a support you can use the characterisation sheet.
2. Once a creamy consistency has been obtained by mixing, we form layers of the denture base material in the JST flask, which has been isolated with ISO-K. Characterisation begins by modelling the area of the marginal gingiva. We then proceed by coloring the marginal gingiva, the attached gingiva and the alveolar musosa.

**Important:** In the areas where the different acrylics meet they must be blended into one another in order to avoid abrupt changes in color.

3. The 8 intensive colors can be mixed in as required on the basis of the individual sketch.

**Important:** Dissolve the intensive colors in the AESTHETIC LT Modelling Monomer first.

4. Fibres (veins) which have been mixed in with the monomer are added after forming the layers of the alveolar mucosa. The required minimum layer thickness is 1 mm. Thin edges which dry out and turn white have to be removed with a scalpel.

---

EN

### 5. Heat-curing in the JST flask

The acrylic dough is placed in the flask and pressed onto the characterised acrylic.

### 6. Cold-curing in the PolyMaster

Pour the mixed AESTHETIC BLUE into the bottom section, allowing it to flow into the gingivobuccal fold and the palate of the model.

7. In order to avoid compressing the characterisations, we close the screw-fit JST flask with a scissor-like movement so that the excess material is pressed out towards the back (dorsal area).

8. The screw opposite the anterior region of the denture is then carefully tightened with the allen key.

### Polymerisation

Please adhere to the instructions for use of AESTHETIC BLUE and AESTHETIC RED acrylics.

### Note

The instructions for use relating to AESTHETIC BLUE and AESTHETIC RED form an integral part of this operatory procedure.

### **Storage instructions**

- Keep out of reach of children
- Note storage instructions and expiry date on the primary packaging
- Do not use the product after the stated expiry date
- Store in an ambient temperature of 2–28 °C / 36–82 °F
- Protect against heat and sunlight

### **Warning**

This material has been developed solely for use in dentistry and must be handled strictly in accordance with the instructions for use. The manufacturer cannot be held liable for any damages arising as a result of failure to observe the instructions for use or the stated area of application. The same applies in the event that the product is mixed or processed with other manufacturers' products. The user shall be solely responsible for testing the material with respect to its suitability prior to use for any purpose other than those explicitly stated in the instructions.

**Safety data sheets can be found on our website at [www.candulor.com](http://www.candulor.com).**

### Danger warnings

- The monomer contains methyl methacrylate (MMA)
- MMA is an irritant and easily flammable (flash point: +10 °C)
- MMA and its vapours are irritating to the eyes, skin and respiratory system
- May cause sensitisation upon skin contact
- Do not inhale vapours
- Keep away from sources of ignition – no smoking

- Do not empty into drains
- Avoid contact of the skin with monomer and uncured material. Many commercial gloves, e.g. those made of latex or vinyl, are not monomer-resistant and therefore do not provide protection against the sensitising effect of methacrylates.
- Wear a mask when grinding and use a suction removal system.
- Always use a spatula when handling the mixture.
- Take precautionary measures against static discharges.

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### AESTHETIC Intensive Opaque / AESTHETIC Intensive Colors

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🏢 Candulor AG