**COLD**

# AUTOPLAST

POLYMER/MONOMER

**INSTRUCTIONS FOR USE**

---

GEBRAUCHSINFORMATION

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---

Dear customer,

The product that you have purchased qualifies as a medical device under European Directive 93/42/EEC. Please ensure that you file the lot number and product name in your administrative system on receipt of the product. For each piece of work that you produce, please make a note of all the materials used along with all the appropriate lot numbers. For your own protection as well as for the protection of your patients and the environment, please also follow the important guidelines in the following instructions for use.

#### **Product description**

AUTOPLAST denture base is a cold-curing, PMMA-based acrylic for the pouring technique. It has excellent flow and modelling properties which make it easy to use. To achieve the best results, please follow the instructions given in this leaflet.

#### **Composition / Powder**

Polymethyl methacrylate, catalyst, pigments

#### **Composition / Liquid**

Methyl methacrylate stab., dimethacrylate, catalyst, stabilizer

### Storage instructions

Store the material in a cool, dark, well-ventilated place.

Storage temperature: 2–28 °C / 36–82 °F. Do not use the material after the expiry date.

Keep out of the reach of children.

### Indication

- Full dentures
- Partial dentures
- Combination dentures
- Dental implants
- Repairs

### Contraindication

Avoid direct contact with unpolymerised material within the oral cavity. If a patient is known to be allergic to any of the ingredients in AUTOPLAST, the material must not be used.

### Side effects

No systemic side effects have been reported to date. In individual cases, local allergic reactions to PMMA/MMA-based denture base materials have been reported.

### Danger warnings

- Monomer contains methyl methacrylate (MMA)
- MMA is an irritant and easily flammable (flash point: + 10 °C / 50 °F)
- MMA and its vapours are irritating to the eyes, skin and respiratory system
- May cause sensitisation upon skin contact
- Do not inhale vapours
- Keep away from sources of ignition – no smoking

- Do not empty into drains
- Avoid contact of the skin with monomer and uncured material. Many commercial gloves, e. g. those made of latex or vinyl, are not monomer-resistant and therefore do not provide protection against the sensitising effect of methacrylates.
- Wear a mask when grinding and use a suction removal system.
- Always use a spatula when handling the mixture.
- Take precautionary measures against static discharges.

## Warning

This material has been developed solely for use in dentistry and must be handled strictly in accordance with the instructions for use. The manufacturer cannot be held liable for any damages arising as a result of failure to observe the instructions for use or the stated area of application. The same applies in the event that the product is mixed or processed with other manufacturers' products. The user shall be solely responsible for testing the material with respect to its suitability prior to use for any purpose other than those explicitly stated in the instructions.

**Safety data sheets can be found on our website at [www.candulor.com](http://www.candulor.com).**

## INSTRUCTIONS FOR USE

### 1. Preparation


After boiling out the water-saturated plaster, isolate the surfaces with two applications of ISO-K separating liquid and allow to dry thoroughly. To ensure proper bonding with the denture base, roughen the teeth well and fit retentions, then wet with monomer.

### 2. Dosage

Ideal mixing ratio for one denture:

13 g polymer : 10 ml monomer

**PolyMaster:** 26 g polymer : 20 ml monomer

The dosage system ensures an ideal mixing ratio and minimum polymerisation shrinkage of AUTOPLAST. The scale on the monomer cylinder is in millilitres, on the polymer beaker in grams.

Too high a monomer content will alter the shade, consistency and handling properties of the material. This may also inhibit polymerisation, possibly leading to irritation of the mucous membranes and to porosity.

### **3. Mixing**

Mix together the powder and liquid, based on the recommended ratio, for approx. 20–30 seconds using the spatula. Blend thoroughly. Leave to stand for 15 seconds so that any air bubbles can escape.

**Mix the powder and liquid thoroughly!**

### **4. Pouring phase**

The pouring phase lasts approximately 2.5–3 minutes at room temperature (18–25 °C / 64–77 °F). Pour the material into the mould within this time span.

## 5. Modelling phase

After a transition phase of roughly 5–6 minutes the material will assume a firmer consistency and can be modelled for a period of 3 minutes. The higher the room temperature, the shorter the pouring and setting phases.

**The synthetic material must be placed in the pressure vessel within no more than 10–11 minutes.**

## 6. Polymerisation

Polymerisation takes place in the pressure vessel for 15 minutes at a temperature of 40 °C and a pressure of 2 bar. The stated polymerisation times must be observed.

If the ideal mixing ratio and the recommended methods are used, a residual monomer content of < 4.5 % can be achieved.

## 7. Finishing

Remove the precast. Check the occlusion of the denture and finish in the usual manner.

### Repairs and corrections of AUTOPLAST

Repairs and corrections may be carried out using either AUTOPLAST itself or AESTHETIC BLUE. The contact surfaces must be well roughened and wetted with monomer. Then apply the mixed material and polymerise as described in the instructions.

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🏢 Candulor AG