Instructions for Use Instrument Set

Technical

- Sliding Caliper
- Papillameter
- Profile Compass
- Rim Former
- Bite Fork
- Modelling Instrument
- Artikont Articulating Ribbon
- ToothScout

Clinical

- Bite Fork Set
- CRS Set 10
- Papillameter
- Sliding Caliper
- Rim Former
- ToothScout

For dental use only

CAUTION: US Federal law restricts this device to sale by or on the order of a licensed dentist.

Papillameter

The Papillameter provides the ideal tool for determining the length of the upper lip and thus the lip closure line. It is also extremely simple to use. The information obtained enables the dental technician to establish the length and position of the anteriors with a considerably higher degree of accuracy. Ask the patient to stand so that you can perform the measurement with the Papillameter. Hold the Papillameter vertically for this purpose. Ask the patient to relax their upper lip. Hold the Papillameter up to the papilla incisiva. The exact length can then simply be read off the instrument. Repeat this process once or twice to ensure a reliable measuring result. Sterilisation in autoclave (134 °C / 273 °F).

Sliding Caliper

The Sliding Caliper is a tool for evaluating the right tooth mould for the patient. Practical experience gained over a long period of time has shown that there is a relationship between the patient's nose width and the shape of the teeth.

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
Clinical

- Set Fourchette occlusale
- Pied à coulisse
- Papillomètre
- Rim Former
- CRS Set 10
- ToothScout

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

Artikont Articulating Ribbon red

This double-sided articulating ribbon is made of silk and marks on both sides. The special dyeing technique applied and its tear-proof design ensure prolonged use and make this ribbon ideal for prosthetic dentistry. Other colours: blue and green. (80 µm)

Modelling Instrument

Double-ended modelling instrument with aluminium handle for preliminary shaping and fine carving. Excellent for shaping the vestibule in wax or for mixing and packing acrylics.

Profile compasses

In-depth analysis of the model is an absolute must for the set-up of full or partial dentures. Tracing the contour of the alveolar ridge is also an important part of accurate model analysis. The simplest way of doing this is to draw the outline on the plaster base. Reproduction of the alveolar ridge ensures the correct placement of the teeth as the contour is still easily identified even after the wax template has been applied. The compasses are used as follows: the curved end of the compasses is set on the alveolar ridge. As this end follows the contour the pencil reproduces a parallel outline on the base of the model. With the aid of the pencil line it is then possible to establish the lowest point of the alveolar ridge, thus allowing exact definition of the masticatory centre.

CRS Set 10

Set for intraoral registration to determine the vertical and centric relation for full dentures. The complete CRS Set 10 is supplied in an aluminium package with a detailed video showing its use.

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THE NAVIGATOR

Natural prosthetics with three-dimensional movement.

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