ULTRA TRAY

LICHTHÄRTENDER KUNSTSTOFF
LIGHT-CURING RESIN
RÉSINE PHOTOPOLYMÉRISABLE
RESINA FOTOINDURENTE
LICHTUITHARDEND LEPELMATERIAAL
RESINA FOTOPOLIMERIZABLE
СВЕТОТВЕРЖДАЕМАЯ ПЛАСТМАССА

DE GEBRAUCHSINFORMATION

EN INSTRUCTION FOR USE

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ULTRA TRAY – LIGHT-CURING RESIN EN

Ultra Tray is a ready-to-use, light-curing resin.

The light-proof container permits optimum utilization of the material. Surplus material can be stored away from light ready for reuse.

Composition

Ultra Tray consists of acrylated oligomers, glass fillers, dimethacrylate, silicon dioxide, catalysts and pigments.

Indications

Ultra Tray can be used for the fabrication of:

- customized trays
- bite plates
- base plates for wax try-ins
- base plates for fixing registration devices
- check bites for framework try-ins

Contraindications

Avoid direct intraoral contact with unpolymerized Ultra Tray material.

Side effects

In patients with a hypersensitivity, repeated contact with unpolymerized material may cause irritation or, in rare cases, lead to a sensitizing reaction. It should not be used for patients with a known allergy to any of the material's ingredients.

Benefits

- quick, clean application
- exact fit
- excellent adaptability
- high stability and good elastic recovery
- uniform layer thickness
- economical, thanks to reuse of surplus material

Instructions for use

After blocking out undercut areas, soak the model in water and isolate with ISO-K. Mould Ultra Tray onto the model using a suitable instrument.

Working time

In normal daylight without direct exposure to sunlight: approx. 10–15 minutes (take the lighting in your working area into consideration).

Adhesion of material to instruments can be prevented with a little water. Cut the functional edge to length with a sharp instrument and remove excess. Use the surplus material to build the tray grip. Residual Ultra Tray material can be reused at a later date. Please keep in mind that exposure to ultraviolet light for an extended period of time can trigger polymerization!

Polymerize the completed tray in a suitable polymerization unit. Any commercially available ultraviolet or halogen light-poly-

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the instructions for use. Liability cannot be accepted for damages resulting from failure to observe the instructions for use or the stipulated area of application. The user is responsible for testing the material with the respect to its suitability prior to use for any purpose not explicitly stated in the instructions. Descriptions and data constitute no warranty of attributes and are not binding. Moreover, liability cannot be accepted if Ultra Tray is mixed with competitors' products or processed in combination with them.

merization unit can be used. Since the light output intensity of such units varies considerably, the material should be cured from both sides without vacuum for 3–5 minutes. Finishing can be done immediately after curing.

Optimum physical properties are achieved after a polymerization time of 5 minutes. Use alcohol to clean the polymerized surface. If necessary, the functional edge can be trimmed with an abrasive instrument. Because Ultra Tray produces a very fine dust when ground it is advisable to work with a suction device and a mask. Use a pumice powder slurry for polishing.

The tray made with Ultra Tray can be used immediately for impression taking.

Modifications and corrections

Roughen or sandblast the polymerized tray/base plate at the area requiring correction with a tungsten carbide bur and clean thoroughly. Then remove new Ultra Tray material from the lightprotected tin, adapt and polymerize again.

Impression materials

All impression materials can be used with the appropriate adhesive – provided that they do not require moist storage.

Suitable polymerization units

- Ivoclar Targis Power – 10 min. program 2
- Any other commercially available ultraviolet or halogen light-polymerization unit

Important

- Place excess material back into the lightprotected container immediately!
- Material demonstrating signs of polymerization should be discarded!
- Remove the soft inhibited surface layer with alcohol!

Storage

- Store protected from light at room temperature
- Storage temperature: 2–28 °C / 36–82 °F
- Do not use the material after the indicated expiry date of expiration
- Keep out of children's reach

Delivery Form

Assortment – Standard pack:
50 maxillary plates (pink) +
1 light-protected container

Risks and dangers

Unpolymerized Ultra Tray material may have an irritating effect. Avoid contact with skin mucous membrane and eyes. Direct skin contact with unpolymerized material may lead to a sensitizing reaction. A mask should be worn and a suction device used during grinding.

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