# CRS SET 10

REGISTRIERSET/REGISTRATION SET

INSTRUCTIONS FOR USE

---

GEBRAUCHSINFORMATION

## INHALTSVERZEICHNIS / CONTENTS

---

		SEITE / PAGE	
DE	Gebrauchsinformation	4	<u>DE</u>
EN	Instructions for use	24	<u>EN</u>
FR	Mode d'emploi	44	<u>FR</u>
IT	Istruzioni d'uso	64	<u>IT</u>
ES	Instrucciones de uso	84	<u>ES</u>
NL	Gebruiksaanwijzing	104	<u>NL</u>
RU	Инструкция по применению	124	<u>RU</u>

## CRS SET 10 REGISTRIERSET

---

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## CRS SET 10 REGISTRATION SET

---

Dear customer,  
We are delighted that you have decided to purchase the CANDULOR CRS Set 10 bite registration set. This registration set is a technically advanced high quality product manufactured according to state-of-the-art technology. However, if used incorrectly, this can lead to risks. Please observe the notes and read the instructions for use.

We hope you enjoy working with the product and wish you every success.

### Notes on the instructions for use

Applicable registration: CRS Set 10

Target group: dentists, dental professionals, dental technicians

These instructions for use serve the safe, appropriate and economical use of the registration.

In case of loss, the instructions for use can be obtained from the service point or downloaded as PDF from the Internet at [www.candulor.com](http://www.candulor.com).

EN

---

## Safety notes

We recommend this section to all persons working with the registration or performing maintenance and cleaning work on the registration.

### Use as intended

When manipulating with the CRS Set 10 on the patient, always ensure that the patient's tongue is not injured by the writing tip of the pen.

### Risks and hazards

#### Aspiration and suffocation hazards

As the registration includes small parts which can become loose if not fixated properly, there is a risk that these parts can be swallowed.

## Injury hazard

To avoid injury due to the pen, the patient must be instructed to place his/her tongue under the pen plate.

### Safety cautions and danger warnings

The CRS Set 10 may only be used exclusively for the scope described in the following section. The intended use also includes:

- > Observing the instructions, regulations and notes in these instructions for use.
- > Correct maintenance and cleaning of the registration.

### **Description of design and functions**

The CRS Set 10 is an intraoral support pin registration device for the individual determination of the centric jaw position in edentulous patients.

The basic principle of the intraoral support pin method is based on McGrane's concept of a three-point support of the two jaw joints and the registration pin, which is placed in the center of the mandible. The support pin draws the typical image of an arrow point on a colored registration plate when the patient performs protrusion or limit movements laterally to the left and right from a dorsal position. The desired central position can then be fixated with the fixa-

tion plate and then the entire registration can be coded three-dimensionally with a suitable material (e.g. registration silicone or plaster). Then the registration can be forwarded to the laboratory for further processing.

### **Indication**

The registration set is used to record the joint-related central position of the mandible for full dentures.

### **Contraindication**

The registration cannot be used for dentulous or partly dentulous situations.


## Scope of delivery

Please check whether the scope of delivery is complete.

### 662513 – CRS Set 10

- > 662514 – Registration tip short ..... 1 piece
- > 662515 – Registration tip long ..... 1 piece
- > 662516 – Registration plate ..... 1 piece
- > 662517 – Pin supporting plate small ..... 1 piece
- > 662518 – Pin supporting plate large ..... 1 piece
- > 662519 – Screw and washer for  
registration plate ..... 1 piece
- > 662520 – Plexiglas plate ..... 1 piece
- > 696276 – CANDULOR Marker ..... 1 piece
- > Instructions for use ..... 1 piece

## CRS SET 10 – COMPONENTS

Please clean the individual components before initial use.

### **Processing the CRS Set 10 registration**

- 1.** The anatomically formed bite templates are prepared in the laboratory and given to the dentist.
- 2.** The back of the Papillameter scale has a small support, which is placed exactly on the incisal papilla of the patient. The measurement is taken using the scale while the upper lip is relaxed. This procedure should be repeated one to two times to achieve an accurate measurement result. This is then used to determine the wax wall height.
- 3.** The Papillameter is placed on the Papilla incisiva of the master model and the measurements taken from the

patient are transferred to dividers or Sliding Calipers. The base edge is best suited as reference point on the model. This measurement is transferred to the wax wall height.

- 4.** The maxilla template fabricated in the laboratory is placed in the mouth and checked to ensure it fits properly. The interpupillary line and Camper's level are checked using a Bite Fork and, if necessary, corrected using the Rim Former. The offset edge of the Rim Former is placed against the tuber tapers. The heated Rim Former is pressed forwards against the wax wall. This melts the entire wax wall parallel to the desired height.

5. Mark one point on the patient's chin and nose tip and measure the resting balance position with the Sliding Caliper. The patient must sit upright and keep his/her straight for bite registration.

6. The mandible template fabricated in the laboratory is placed in the mouth and checked to ensure it fits properly. Here, the wax wall height should match up with the equator of the tongue.

7. Check the position of the wax wall templates with respect to each other as well as the vertical height. It should be 2 – 3 mm lower than the relaxed rest position. Drawing the center line on the maxilla.

8. Using these measurements, it is recommend making a temporary articulation of the model.

9. The mandible registration plate is fixed to the occlusal level and notches are made in the wax wall to aid keying (Fig. 1).


Fig. 1

**10.** The supporting pin plate is heated and immersed in the wax wall until the registration pin in regio 15/25 is positioned at the height of the occlusal level and lies vertical. The wax wall is to be reduced such that there is no longer any contact between the upper and lower wax wall during the movements. Create notches to help with keying (Fig. 2).

Fig. 2

**11.** Place the registration template in the mouth of the patient. Only the registration pin and registration template may touch each other during the registration procedure. While sitting in the upright position, the patient is instructed to move his mandible as far to the left as possible and then back to the center while the pin and registration plate are touching. This procedure is also repeated on the right side two to three times. The patient is then asked to move his mandible to the front and back. This exercise is to prepare the patient for the actual bite registration procedure.

**Attention:** to avoid injury due to the pen, the patient must be instructed to place his/her tongue under the pen plate or in dorsal direction.

**12.** Using the CANDULOR Marker provided, generously paint the mandibular registration plate. The templates are then put in and, following the dentist's instructions, the patient traces the angle by anterior, posterior and lateral movement of the mandible (Fig. 3).

**13.** Crosshair symbol drawn on the tip of the arrow (centric).


Fig. 3

**14.** To check the gothic arch, the entire plate except for the crosshairs is painted again and the entire bite registration procedure is repeated. If the second gothic arch is also positioned exactly on the crosshairs, the Plexiglas plate can be placed in the position of the tip of the arrow fixed into place using the screw and washer (Fig. 4).


Fig. 4

**15.** The opening and closing movements make it possible to check the mouth and see whether the registration pin snaps into the hole of the Plexiglas plate.

**16.** If this is the case, the position is cemented using registration silicone or impression plaster.

**17.** If this is not the case, the registration procedure must be repeated.

### Notes on cleaning

What	When	With what
Registration parts and screws (metal)	After each use	Sterilization in the autoclave at 134° C for 5:30 minutes, perform pre-vacuum method. Please observe the applicable standards.
Plexiglas plate (acrylic)	After each use	Perform disinfection in the immersion bath according to the instructions of the disinfectant manufacturer.

### Notes on care

Avoid any contact with strong acids and solvents (e.g. MMA, acetone). This prevents any damage to the surface treatment.

Clean, sterilize or disinfect the parts after every use or contact with saliva or blood.

- > Only use disinfectants suitable for stainless steel and light metals.
- > Avoid any contact with strong acids, alkalis or solvents.
- > Clean with ultrasonic bath, water, steam jet or sterilization equipment.

### Warning

The registration has been developed for dental use. Initial operation and handling must be conducted according to the instructions for use. CANDULOR does not accept liability for damages caused by use other than directed or improper handling. The user is also responsible for confirming that the registration is suitable for the intended purpose before use, particularly if the purpose is not listed in the instructions for use.

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**🏢 CANDULOR AG**