



AESTHETIC RED / AESTHETIC BLUE

Certificate of Biocompatibility

Test material: AESTHETIC RED AND AESTHETIC BLUE Denture Base Material

Manufacturer CANDULOR AG, Boulevard Lilienthal 8, CH-8152 Glattpark (Opfikon), Switzerland

Standard composition
AESTHETIC RED Heat-curing denture basematerial Polymethyl methacrylate with a cross-linked monomer matrix
AESTHEIC BLUE Auto-curing denture basematerial Polymethyl methacrylate with a cross-linked monomer matrix

Classification
AESTHETIC RED denture base materials comply with the ISO 20795-1:2013 standard and, given their chemical composition and polymerization temperature, belong to the heat-curing polymers (Type 1, Class 1).
AESTHETIC BLUE denture base materials comply with the ISO 20795-1:2013 standard and, given their chemical composition and polymerization temperature, belong to the auto-curing polymers (Type 2, Class 1).

Residual monomer content testing according to ISO 20795-1:2013
Limit for auto-curing polymers 4.5 percent
Limit for heat-curing polymers 2.2 percent
Value achieved by AESTHETIC RED ≤ 2.2 percent
Value achieved by AESTHETIC BLUE ≤ 4.5 percent

Water solubility testing according to ISO 20795-1:2013
Limit for auto-curing polymers $\mu\text{g}/\text{mm}^3 < 8.0$
Limit for heat-curing polymers $\mu\text{g}/\text{mm}^3 < 1.6$
Example value for AESTHETIC RED $0.1\mu\text{g}/\text{mm}^3$
Example value for AESTHETIC BLUE $2.7\mu\text{g}/\text{mm}^3$

Cytotoxicity Results: Polymerized denture base material was extracted with culture media for 24 h at 37 °C and the resulting elute was tested in an in vitro cytotoxicity test (XTT). No cytotoxicity could be detected.

Genotoxicity Results: Polymerized denture base material was extracted and the resulting eluate was tested. In both test systems no genotoxicity could be detected.

Glattpark, 19. April 2016

Dr. Mathrin Fischer
Scientific Service

Manufacturer
CANDULOR AG
Boulevard Lilienthal 8
CH-8152 Glattpark (Opfikon)
Switzerland
Tel. +41 (0)44 805 90 00
Fax +41 (0)44 805 90 90
www.candulor.com

Notified Body
TÜV Süd Product Service GmbH
Riedlerstrasse 65
D-80339 München
Germany

CE 0123