



Compatibility of products with disinfectants

29.05.2020 / R&D / Candulor AG

Overview of the compatibility of disinfectants with the dental products and equipment of Candulor.

Material / product	Compatibility with disinfectants				Notes
	70% (v/v) i-propanol (Reaction time: 1 min) ¹	70% (v/v) ethanol (Reaction time: 1 min) ²	0.1% sodium hypochlorite (Reaction time: 1 min) ²	0.5% hydrogen peroxide (Reaction time: 1 min) ²	
Removable partial / complete denture prosthetics - Teeth - Denture base material	✓ suitable	✓ suitable	✓ suitable	✓ suitable	n.a.
Auxiliary materials - Waxes	✓ suitable	✓ suitable	✓ suitable	✓ suitable	n.a.
Tray material - C-Plast	✓ suitable	✓ suitable	✓ suitable	✓ suitable	n.a.
Instruments / registration - CRS Sets - Bite fork	✓ suitable	✓ suitable	✓ suitable	✓ suitable	<i>Disinfection according to instructions for use, it can be autoclaved afterwards</i>

¹ https://vah-online.de/files/download/vah-mitteilungen/2020_19Mrz_VAH%20Mitteilung%20Händedesinfektion%20Rezepturen%20aktualisiert.pdf

² G. Kampf et al. / Journal of Hospital Infection 104 (2020) 246-251